1 Yichao

For shutting off, go through the checklist in the following order

- Turn off equipments
- Turn off power supplies
- Turn off computers

Wait before everyone else finish their corresponding steps before starting the next one.

For powering up, go through the checklist in the following order

- Turn on computers
- Turn on power supplies
- Turn on equipments
- Check computer connection to devices

Also wait before everyone else finish their corresponding steps before starting the next one.

- 2 Equipments
- 3 Power supplies
- 4 Computers

5 Jon

For shutting off, go through the checklist in the following order

- Turn off equipments
- Turn off power supplies
- Turn off computers

Wait before everyone else finish their corresponding steps before starting the next one.

For powering up, go through the checklist in the following order

- Turn on computers
- Turn on power supplies
- Turn on equipments
- Check computer connection to devices

Also wait before everyone else finish their corresponding steps before starting the next one.

- 6 Equipments
- 7 Power supplies
- 8 Computers

9 Lee

For shutting off, go through the checklist in the following order

- Turn off equipments
- Turn off power supplies
- Turn off computers

Wait before everyone else finish their corresponding steps before starting the next one.

For powering up, go through the checklist in the following order

- Turn on computers
- Turn on power supplies
- Turn on equipments
- Check computer connection to devices

Also wait before everyone else finish their corresponding steps before starting the next one.

10 Equipments

11 Power supplies

12 Computers

13 Jessie

For shutting off, go through the checklist in the following order

- Turn off equipments
- Turn off power supplies
- Turn off computers

Wait before everyone else finish their corresponding steps before starting the next one.

For powering up, go through the checklist in the following order

- Turn on computers
- Turn on power supplies
- Turn on equipments
- Check computer connection to devices

Also wait before everyone else finish their corresponding steps before starting the next one.

- 14 Equipments
- 15 Power supplies
- 16 Computers

17 Frederic

For shutting off, go through the checklist in the following order

- Turn off equipments
- Turn off power supplies
- Turn off computers

Wait before everyone else finish their corresponding steps before starting the next one.

For powering up, go through the checklist in the following order

- Turn on computers
- Turn on power supplies
- Turn on equipments
- Check computer connection to devices

Also wait before everyone else finish their corresponding steps before starting the next one.

18 Equipments

19 Power supplies

20 Computers

• Of Andor Camera for 1.0

Turn off check list

- Switch off
- Switch off chiller
- Unplug camera power (zip tie to camera)

Turn on check list

- Plug in camera power
- Switch on chiller
- Switch on

• Off Andor Camera for 1.5 Turn off check list

- I di ii on check in
 - Switch offSwitch off chiller
 - Unplug camera power (zip tie to camera)

Turn on check list

- Plug in camera power
- Switch on chiller
- Switch on

• Of Tweezer 1 computer Check following devices

- Thorlabs cameras
- Na Raman amplifier

• Off Tweezer 2 computer Check following devices

– Andor camera

• Off On Tweezer 3 computer Check following devices

- Software radio

• Tweezer 4 computer Check following devices

- Andor camera
- Thorlabs camera

• PLL for 1.0 Cs MOT

Turn off check list

- Unplug all signals (zip tie together)
- Unplug 15V and 5V powers (zip tie together)

Turn on check list

- Plugin 15V and 5V powers
- Plugin all signals
- Apply settings according to note on Generic tab, Equipment Settings.

• Off PLL for 1.0 Cs Raman

Turn off check list

- Unplug all signals (zip tie together)
- Unplug 15V and 5V powers (zip tie together)

Turn on check list

- Plugin 15V and 5V powers
- Plugin all signals
- Apply settings according to note on Generic tab, Equipment Settings.

• Off On PLL for 1.5 Cs MOT

Turn off check list

- Unplug all signals (zip tie together)
- Unplug 15V and 5V powers (zip tie together)

Turn on check list

- Plugin 15V and 5V powers
- Plugin all signals
- Apply settings ???

• Off On 1.5 Feshbach coils

Before power shutdown

- Wrap exposed part in napkins to prevent condensations

• Off On 1.5 IGBT

Before power shutdown

- Wrap exposed part in napkins to prevent condensations

• Computer control box for 1.0

Turn off check list

- (Do following three steps quickly)
- Turn off the box by flipping the switch on the front side of the box
- Unplug the high voltage (48V) power supply (circular plug) on the front side of the box. Then unplug the power supply itself fromm the outlet.
- Unplug the 12V board power supply. (L.T.E. one with a rectangular connector in the front)
- Unplug USB power on the front side near the fan.
- Zip tie the three unpluged powers to the computer control box.

- Connect to 3.5G Windfrek clock generator from Tweezer1. Set frequency to 3.5G and amplitude to max.
- (Remove power connection zip ties and) Plug in USB power supply.
- − (Do following three steps quickly)
- Plug in the 12V board power supply.
- Plug in high voltage power supply.
- Turn on the power switch of the box. The fans should start spinning.

• Off Computer control box for 1.5 Turn off check list

- (Do following three steps quickly)
- Turn off the box by flipping the switch on the front side of the box
- Unplug the high voltage (48V) power supply (circular plug) on the left side of the box. Then unplug the power supply itself fromm the outlet.
- Unplug the 12V board power supply. (L.T.E. one with a rectangular connector in the front)
- Unplug USB power on the front side near the fan.
- Zip tie the three unpluged powers to the computer control box.

Turn on check list

- Connect to 3.5G Windfrek clock generator from Tweezer1. Set frequency to 3.5G and amplitude to max.
- (Remove power connection zip ties and) Plug in USB power supply.
- (Do following three steps quickly)
- Plug in the 12V board power supply.
- Plug in high voltage power supply. The fans should start spinning.
- Turn on the power switch of the box.

• Off On Homebuilt temperature servo for Cs 1.0 MOT

1st level above Cs tweezer/MOT lasers Switch off/on the front switch

• Off On Benchtop power supply 1st level above Cs tweezer/MOT lasers Record voltage/current values:

Left Voltage	
Left Current	
Right Voltage	
Right Current	

Turn off check list

- Unplug loads
- Label and zip tie loads to the power supply
- Turn off the power supply
- Unplug from wall

Turn on check list

- Plug into the wall
- Turn on the power supply Do NOT do this with loads pluged in
- Plug the load back in

• Ist version Till's current controller for Cs MOT

1st level above Cs tweezer/MOT lasers Record current value:

Current	
---------	--

Turn off check list

- Turn off
- Unplug power (zip tie to controller)

- Plug in power
- Turn on

• Of On Homebuilt $\pm 15,5$ V power supply

2nd level above Cs tweezer/MOT lasers Turn off check list

- Unplug load
- Label and zip tie loads to the power supply
- Switch off on the back

Turn on check list

- Switch on
- Plug the load back in
- Off On Cs 1.0 MOT lock box 1st level above Cs MOT lasers
 - Make sure the laser is unlocked.
- Off Novatech function generator
 1st level above Cs MOT lasers Turn off check list
 - Save settings

Turn on check list

- Restore settings

• New Focus Cs Raman lasers drivers (2x)

1st level above Cs Raman lasers

Record Temperature, Current, Piezo voltage Temperature is accessible in the menu under system status

F3 Current	
F3 Piezo Voltage	
F3 Temperature	
F4 Current	
F4 Piezo Voltage	
F4 Temperature	

Turn off check list

- Turn off
- Unplug power

- Plug in power
- (Turn on)

• Benchtop power supply
1st level next to Cs Raman driver
Record voltage/current values:

Voltage	
Current	

Turn off check list

- Unplug loads
- Label and zip tie loads to the power supply
- Turn off the power supply
- Unplug from wall

Turn on check list

- Plug into the wall
- Turn on the power supply **Do NOT do this with** loads pluged in
- Plug the load back in
- Off On Homebuilt 24V power supply 2nd level above Cs Raman lasers
 Turn off check list
 - Unplug load
 - Label and zip tie loads to the power supply
 - Switch off on the back

Turn on check list

- Switch on
- Plug the load back in
- Off On SRS shutter driver for Cs
 1st level above Cs Raman lasers

Switch off/on with the switch on the back.

• Off New Focus lock box for Cs Raman

1st level above Cs Raman lasers Unplug power.

• Benchtop power supply for Cs Raman locking

2nd level above Cs Raman lockbox Record voltage/current values:

Voltage	
Current	

Turn off check list

- Unplug loads
- Label and zip tie loads to the power supply
- Turn off the power supply
- Unplug from wall

- Plug into the wall
- Turn on the power supply Do NOT do this with loads pluged in
- Plug the load back in

• Of Timebase driver for Na Raman and Na D2 (x2)

1st level above Na Raman and D2 seed lasers Record current and T1 settings

Raman Current	
Raman T1	
MOT Current	
MOT T1	

Turn off check list

- Turn off
- Unplug power

Turn on check list

- Plug in power
- Turn on
- Set current/temperature.
- Enable T1 survo
- Turn on current

• Off On Thorlabs temperature controller for Na Raman doubler

1st level above Na Raman and D2 seed lasers Record temperature settings

Setpoint

Turn off check list

- Turn off
- Unplug power

Turn on check list

- Plug in power
- Check temperature setting
- Turn on

• of Daniebuilt lockbox for Na D2

1st level above Na seed

Turn on-off together with the adjustable high voltage power supply above

Record high voltage power supply output

Voltage	
---------	--

Turn off check list

- Unlock
- Turn down the voltage of the high voltage power supply to 0
- Unplug high voltage power on lockbox (zip tie)
- Unplug 15V from lockbox (zip tie)
- Turn off high voltage power supply
- Unplug high voltage power supply from wall

Turn on check list

- Plug high voltage power supply into the wall
- Turn on high voltage power supply (make sure the output is 0)
- Plug in 15V to lockbox
- Plug in high voltage to lockbox
- Turn up the voltage of the high voltage power supply
- Check HV output

• Oscilloscope for Na D2 lock 1st level above Na Raman amplifier

– Unplug from wall

•	Off	On	MPB	Ra Ra	aman	ampl	lifier
						am path	
	Turn	off	check li	\mathbf{st}			

- Power off
- Unplug from wall

Turn on check list

- Plug into wall
- Power on

• Off On Variac for Na cell Above Na D2 beam path

Record temperature and voltage settings

Voltage	
Temperature	

Turn off check list

- Turn voltage to 0
- Switch off
- Unplug from wall

Turn on check list

- Plug into wall
- Switch on
- Turn the voltage up **SLOWLY** and monitor the temperature at the same time.

• Off On SRS shutter driver for Na 1st level above Na D2 beam path, below Variac Switch off/on with the switch on the back.

• Off Denchtop power supply

2nd level above Na Variac

Record voltage/current values:

Left Voltage	
Left Current	
Right Voltage	
Right Current	

Turn off check list

- Unplug loads
- Label and zip tie loads to the power supply
- Turn off the power supply
- Unplug from wall

- Plug into the wall
- Turn on the power supply Do NOT do this with loads pluged in
- Plug the load back in

• Off On Homebuilt lockbox for Na D1

1st level above Na D2 beam path Turn on-off together with the fixed high voltage power supply above

Turn off check list

- Unlock
- Unplug high voltage power on lockbox (zip tie)
- Unplug 15V from lockbox (zip tie)
- Unplug high voltage power supply from walls

Turn on check list

- Plug high voltage power supply into the wall
- Plug in 15V to lockbox
- Plug in high voltage to lockbox
- Check HV output

• Off On Thorlabs temperature controller for Na D1 doubler

1st level above Na D1 Record temperature settings

Setpoint	

Turn off check list

- Turn off
- Unplug power

Turn on check list

- Plug in power
- Check temperature setting
- Turn on

• of Thorlabs temperature controller for Na D1 seed

1st level above Na D1 seed Record temperature settings

Setpoint	

Turn off check list

- Turn off
- Unplug power

Turn on check list

- Plug in power
- Check temperature setting
- Turn on

• Of Thorlabs current controller for Na D1 seed

1st level above Na D1 seed

Record current settings

Current

Turn off check list

- Turn off
- Unplug power

- Plug in power
- Check current setting
- Turn on

• Off Denchtop power supply 2nd level above Na D1 seed controllers
Record voltage/current values:

Left Voltage	
Left Current	
Right Voltage	
Right Current	

Turn off check list

- Unplug loads
- Label and zip tie loads to the power supply
- Turn off the power supply
- Unplug from wall

Turn on check list

- Plug into the wall
- Turn on the power supply Do NOT do this with loads pluged in
- Plug the load back in
- Homebuilt 24V power supply 2nd level above Cs tweezer controllers
 Turn off check list
 - Unplug load
 - Label and zip tie loads to the power supply
 - Switch off on the back

Turn on check list

- Switch on
- Plug the load back in

• Off On Benchtop power supply
2nd level above Cs tweezer, right of 24V power supply

Record voltage/current values:

Left Voltage	
Left Current	
Right Voltage	
Right Current	

Turn off check list

- Unplug loads
- Label and zip tie loads to the power supply
- Turn off the power supply
- Unplug from wall

Turn on check list

- Plug into the wall
- Turn on the power supply Do NOT do this with loads pluged in
- Plug the load back in
- In Homebuilt temperature servo for Cs tweezer

1st level above Cs tweezer Switch off/on the front switch

•	Off	On	2nd	version	Till's	current
				or Cs Tw		

1st level above Cs tweezer

Record current value:

Current

Turn off check list

- Turn off
- Unplug power (zip tie to controller)

Turn on check list

- Plug in power
- Turn on

• Of Thorlabs current controller

1st level below Cs tweezer driver Unused

Turn off check list

- Turn off
- Unplug power

• Of Thorlabs temperature controller for Cs tweezer

1st level above Cs tweezer Record temperature settings

Setpoint	

Turn off check list

- Turn off
- Unplug power

Turn on check list

- Plug in power
- Check temperature setting
- Turn on

• Off Intensity servo for Cs
1st level above Cs tweezer
Unplug power

• Off On Benchtop power supply

1st level above Cs tweezer

Record voltage/current values:

Left Voltage	
Left Current	
Right Voltage	
Right Current	

Turn off check list

- Unplug loads
- Label and zip tie loads to the power supply
- Turn off the power supply
- Unplug from wall

- Plug into the wall
- Turn on the power supply Do NOT do this with loads pluged in
- Plug the load back in

•	Off On	2nd	version or Cs 1.5	Till's	current
	COHUIC	TICI I	$\mathbf{o}_{\mathbf{i}}$	(A = I)	
	1st level	above 1	$1.5~\mathrm{Cs}~\mathrm{MOT}$, ,	

Record current values:

Top Current	
Bottom Current	

Turn off check list

- Turn off
- Unplug power (zip tie to controller)

Turn on check list

- Plug in power
- Turn on
- Off On Thorlabs current controller
 1st level below Cs 1.5 MOT drivers
 Unused
 Turn off check list
 - Turn off
 - Unplug power

• Denchtop power supply 2nd level above Cs 1.5 MOT controllers

Record voltage/current values:

Left Voltage	
Left Current	
Right Voltage	
Right Current	

Turn off check list

- Unplug loads
- Label and zip tie loads to the power supply
- Turn off the power supply
- Unplug from wall

Turn on check list

- Plug into the wall
- Turn on the power supply Do NOT do this with loads pluged in
- Plug the load back in
- In the Homebuilt 24V power supply 2nd level above Cs 1.5 MOT controllers

 Turn off check list
 - Unplug load
 - Label and zip tie loads to the power supply
 - Switch off on the back

- Switch on
- Plug the load back in

• Off On Homebuilt temperature servo for Cs 1.5 MOT (x2) 1st level above Cs 1.5 MOT

1st level above Cs 1.5 MO: Switch off/on the front switch

• Off On Oscilloscope for Cs 1.5 MOT lock

1st level above Cs 1.5 MOT

- Unplug from wall
- $\stackrel{\text{Off}}{\longrightarrow}$ Homebuilt lockbox for Cs 1.5 $\stackrel{\text{MOT}}{\longrightarrow}$

1st level above Cs 1.5 MOT Turn off check list

- Unlock
- Unplug 15V from lockbox (zip tie)

Turn on check list

- Plug in 15V to lockbox
- Of Function generator for Cs MOT 1.5 lockbox

1st level under C
s MOT 1.5 lockbox Turn off check list

- Record settings
- Turn off
- Unplug from wall

Turn on check list

- Plug into wall
- Turn on
- Restore settings

• Off Denchtop power supply (x2)
1st level above Na 1.5 bean path

Record voltage/current values:

Left Voltage	
Left Current	
Right Voltage	
Right Current	

Turn off check list

- Unplug loads
- Label and zip tie loads to the power supply
- Turn off the power supply
- Unplug from wall

Turn on check list

- Plug into the wall
- Turn on the power supply Do NOT do this with loads pluged in
- Plug the load back in
- If the interval $\pm 15,5$ We have supply $\pm 15,5$ We have $\pm 15,5$ We ha

2nd level above Na 1.5 beam path Turn off check list

- Unplug load
- Label and zip tie loads to the power supply
- Switch off on the back

- Switch on
- Plug the load back in

• Of the HV amplifiers for 1.0 piezo mirrors

Above the space beween Na 1.5 beam path and NaCs 1.5 chamber

Turn on-off together with the function generator below and the high voltage power supply on the left

Turn off check list

- Turn off the function generator
- Unplug function generator output (zip tie and label)
- Unplug HV power supply from the wall

Turn on check list

- Plugin HV power supply to the wall
- Plugin function generator output
- Turn on the function generator

• Off Description Benchtop power supply (x2)

2nd level above piezo mirror electronics Record voltage/current values:

Left Voltage	
Left Current	
Right Voltage	
Right Current	

Turn off check list

- Unplug loads
- Label and zip tie loads to the power supply
- Turn off the power supply
- Unplug from wall

Turn on check list

- Plug into the wall
- Turn on the power supply Do NOT do this with loads pluged in
- Plug the load back in

• Off On Pump

1st level above NaCs 1.0 chamber Turn off check list

- Turn off
- Unplug from walls

- Plug into wall
- Turn on (should be enough)

• Off Denchtop power supply 2nd level above NaCs 1.0 chamber Record voltage/current values:

Left Voltage	
Left Current	
Right Voltage	
Right Current	

Turn off check list

- Unplug loads
- Label and zip tie loads to the power supply
- Turn off the power supply
- Unplug from wall

Turn on check list

- Plug into the wall
- Turn on the power supply Do NOT do this with loads pluged in
- Plug the load back in
- Off On Homebuilt 24V power supply 2nd level above NaCs 1.0 chamber, facing the drawers/wall

Turn off check list

- Unplug load
- Label and zip tie loads to the power supply
- Switch off on the back

Turn on check list

- Switch on
- Plug the load back in

• Off On Benchtop power supply for dispensers

1st level above NaCs 1.0 chamber Turn off check list

- Unplug loads
- Label and zip tie loads to the power supply
- Turn off the power supply
- Unplug from wall

Turn on check list

- Plug into the wall
- Turn on the power supply Do NOT do this with loads pluged in
- Plug the load back in

🔸 🍱 💴 Quantel Laser

1st level in the middle, facing 1.5 table Ask Frederic.

• Off On Intensity servo for Na
1st level above Na tweezer, below 1.0 computer control box
Unplug power

• Off On Benchtop power supply for 1.5 computer control clock

2nd level next to 1.5 computer control box Record voltage/current values:

Voltage	
Current	

Turn off check list

- Unplug loads
- Label and zip tie loads to the power supply
- Turn off the power supply
- Unplug from wall

Turn on check list

- Plug into the wall
- Turn on the power supply Do NOT do this with loads pluged in
- Plug the load back in
- Off on 1st version Till's current controller for butterfly laser

1st level above TiSapph Turn off check list

- Turn off
- Unplug power (zip tie to controller)
- Of Homebuilt temperature servo for butterfly laser

1st level above TiSapph Switch off/on the front switch

• Of On Homebuilt $\pm 15,5$ V power supply

2nd level above TiSapph Turn off check list

- Unplug load
- Label and zip tie loads to the power supply
- Switch off on the back

Turn on check list

- Switch on
- Plug the load back in

• Off On TiSapph

Turn off check list

- Turn off pump, unplug from wall
- Turn off controller, unplug from wall
- Turn off chiller, unplug from wall

- Plug into wall, turn on chiller
- Plug into wall, turn on controller
- Plug into wall, turn on pump

• Off On Agilent function generator for Na 1.0 switching

1st level above TiSapph facing computers Record parameters:

Ch1 Hi	
Ch1 Lo	
Ch2 Hi	
Ch2 Lo	

Turn off check list

- Turn off
- Unplug from wall

Turn on check list

- Plug into wall
- Turn on
- Set frequency, duty cycle, high/low voltages
- Set switching phase back

• Of the HV amplifiers for 1.5 piezo mirrors

Top level of 1.5 rack

Turn on-off together with the function generator below and the high voltage power supply above the 1.0 machine table

Turn off check list

- Turn off the function generator
- Unplug function generator output (zip tie and label)
- Unplug HV power supply from the wall

Turn on check list

- Plugin HV power supply to the wall
- Plugin function generator output
- Turn on the function generator

• Off Denchtop power supply

2nd level (from the top) 1.5 Rack. Next to coil servo

Record voltage/current values:

Left Voltage	
Left Current	
Right Voltage	
Right Current	

Turn off check list

- Unplug loads
- Label and zip tie loads to the power supply
- Turn off the power supply
- Unplug from wall

- Plug into the wall
- Turn on the power supply Do NOT do this with loads pluged in
- Plug the load back in

• Off On Benchtop power supply
Bottom level 1.5 Rack. Next to pump
Record voltage/current values:

Voltage	
Current	

Turn off check list

- Unplug loads
- Label and zip tie loads to the power supply
- Turn off the power supply
- Unplug from wall

Turn on check list

- Plug into the wall
- Turn on the power supply Do NOT do this with loads pluged in
- Plug the load back in

• Off On Pump

Bottom level 1.5 Rack Turn off check list

- Turn off
- Unplug from walls

Turn on check list

- Plug into wall
- Turn on (should be enough)

• Off Cooling water

Valves to close/reopen: